

SPECIFICATION SHEET

Product Name: G-17-Ag-BSA **Catalog #:** L5B01001

Description: BSA Conjugated Human GASTRIN 17 (G-17)

Format: liquid, G-17 purified by HPLC

Purification: Purified

Concentration& Lot.: 10mg/ml **Lot:** L190901

Buffer: 10Mm PBS, pH 7.5

Preservative: None

Description: Endogenous peptide, mainly found in the gut, that activates gastric acid secretion, delays gastric emptying and increases blood flow in the gastric mucosa.

Storage&Shelf life: Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.