

**SPECIFICATION SHEET**

<b>Product Name:</b>	Human CA24-2 Calib. Grade	<b>Catalog #:</b>	L2C00602
<b>Description:</b>	Human Cancer Antigen 24-2 (CA24-2) (Gastrointestinal Cancer) Calibration Grade		
<b>Source:</b>	Human colion carcinoma.		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	Gel filtration and column chromatography.		
<b>Concentration:</b>	450KU/ml (Canag ELISA)	<b>Lot #:</b>	L170201
<b>Buffer:</b>	PBS, pH 7.4		
<b>Preservative:</b>	0.05% Sodium Azide		
<b>Applications:</b>	Immunogen, calibrator in immunoassays and in vitro studies. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage:</b>	Upon receipt, store at -20°C.		
<b>Inactivation:</b>	Not applicable.		
<b>Warnings:</b>	All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		