

SPECIFICATION SHEET

Product Name:	CEA Antigen Grade	Catalog #:	L2C01001
Description:	Carcinoembryonic Antigen (CEA) Antigen Grade		
Source:	Human colon cancer		
Format:	Liquid		
Purification:	Affinity purified (SDS-PAGE 90%)		
Concentration:	1 mg/ml	Lot #:	L180601
Buffer:	0.01M PBS, pH 7.4		
Preservative:	0.05% NaN ₃		
Applications:	Suitable as an immunogen for antisera production, tracer for iodination and standard for immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage:	Short term store at 2~8°C, long term store at -20°C.		
Inactivation:	Not applicable.		
Warnings:	This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.