

**SPECIFICATION SHEET**

**Product name** Trypsinogen-2 (TG2) GST-tag

**Protein description** Recombinant full length Human Trypsinogen-2 (amino acids 1-247) with N terminal GST tag;

**Molecular weight** 54.5 KD

**Protein length** 487 aa

**Expression host** E. coli **Purity** 80% by SDS-PAGE

**Purification notes** Sepharose derived purification.

**Biological activity** Immunoreactive with sera of acute pancreatitis, and cystic fibrosis.

**Form** Liquid

**Storage instructions** Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

**Storage buffer**

**Preservative** None

**Application** ELISA ; WESTERN