

SPECIFICATION SHEET

Description:	ProGRP31-98 antigen	Catalog #:	L2C03002
Source:	Recombinant, E.coli		
MW:	32kd		
Format:	Purified, Lyophilized		
Purification:	Purified, > 90% (determined by SDS gelelectrophoresis)		
Lot:	L161101		
Net Weight:	0.1mg		
Buffer:	Lyophilized from 10mM Tris-HCL, PH8.0 ± 0.2		
Preservative:	None		
Storage:	Store at -20°C. Avoid multiple freeze/thaw cycles.		