

SPECIFICATION SHEET

Product Name: Recombinant 2019-nCov Spike-RBD Protein (2019-nCov-S1-RBD)

Catalog #: L10H002

SYNONYMS: 2019-nCov-Spike-RBD, COVID-19-Spike-RBD,

Source: E. Coli **Species:** Novel coronavirus pneumonia

Molecular weight: 26.5KD (with an 6*his tag at the C-terminus)

Concentration: **Lot #:**

Purity : >95% by SDS-PAGE

Formulation : The Recombinant 2019-nCov-S-RBD is supplied in a 0.2um filtered concentrated solution contains PBS (pH 7.4) .

Application: ELISA; Positive Control; Immunogen; WB;

Storage & Expiration: Upon receipt, aliquot and store at <-20°C .for long term.
Avoid repeated freeze and thaw cycles.

DESCRIPTION: The recombinant 2019-nCoV Spike-RBD Protein consists of the full-length protein 232 amino acids with a 6*his tag at the C-terminus and predicts a molecular mass of 26.5 kDa.