

SPECIFICATION SHEET

Product Name: HIV-M(gp41+gp36+O) **Catalog #:** L2H00101

Description: **HIVgp41+gp36+O Recombinant antigen,**
Antigen Grade

Source: Recombinant ,E.col

Format: Purified, Protein solution in Tris-HCl pH10.0,

Purification: Purity >90%

Concentration&Lot: 5mg/ml **Lot :** L180901

Buffer: 50mM Tris, pH 8.0, 8M Urea, 300mM Nacl

Preservative: None

Application: Immunoassays or Immunization.

Storage&Shelf life: Store at 2~80 C for 1 years. Avoid multiple freeze/thaw cycles.

Warnings: For research only