

SPECIFICATION SHEET

Product Name:	Rabbit anti HIV-I p24	Catalog #:	L2H00501
Description:	Rabbit Antibody to Human Immunodeficiency Virus Type I (HIV-I), p24		
Host Animal:	Rabbit		
Specificity:	HIV-1 p24 (Recombinant antigen)		
Format:	Purified, Liquid		
Concentration&Lot:	2mg/ml	Lot:	L180301
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA, IFA and Western blot. Also suitable for conjugation purposes. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.