

SPECIFICATION SHEET**Product Name:** Recombinant HAdv3 **Catalog #:** L6H001**Description:** Recombinant human Adenovirus Antigen, type III**Background:** Acute diarrhea disease in young children is a major cause of morbidity worldwide and is a leading cause of mortality in developing countries. Adenovirus is the most common agent responsible for acute gastroenteritis, mainly in young children. Adenovirus can infect respiratory tract, gastrointestinal tract, urinary tract, eye and liver. Adenovirus serotype III is known about 1/3 is typically associated with human disease.**Mw:** About 35KD**Source:** Recombinant, E. coli**Format:** Purified, Protein solution in PBS**Purification:** Purity >90%**Concentration&Lot:** 1mg/ml **Lot:** L170701**Buffer:** 10Mm PBS, pH 7.4± 0.02,**Preservative:** None**Application:** Immunoassays as EIA, IFA or Immunization.**Storage&Shelf life:** Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.**Warnings:** For research only