

SPECIFICATION SHEET

Product Name :	Goat Anti HBsAg PcAb	Catalog #:	L1H00201				
Description:	Goat Antibody to Human Hepatitis B Surface Antigen (HBsAg)						
Specificity:	Reacts with both ad and ay antigens. No reaction when tested against normal human serum and single arc when tested against HBsAg elevated human serum by IEP.						
Host Animal:	Goat						
Immunogen:	Purified Human Hepatitis B Surface Ag						
Format:	Purified, Liquid						
Concentration:	3.8 mg/ml	Lot. :	L180601				
Buffer:	10Mm PBS, pH 7.4						
Preservative:	None						
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1H00202</td><td>L1H00201</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1H00202	L1H00201
<u>Capture</u>	<u>Detection</u>						
L1H00202	L1H00201						

Storage& Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.