

SPECIFICATION SHEET

Product Name:	McAb to Myoglobin	Catalog #:	L3C00301				
Description:	Monoclonal Antibody to Human Myoglobin						
Specificity:	Myoglobin						
Clone:	MY001						
Host Animal:	Mouse						
Isotype:	IgG						
Source:	Ascites						
Immunogen:	Recombinant Human Myoglobin						
Format:	Purified, Liquid						
Purification:	Protein G chromatography						
Concentration&Lot:	2.7 mg/ml	Lot:	L190601				
Buffer:	10mM PBS, pH 7.4						
Preservative:	none						
Applications:	Can be used in an immunoassay for specific determination of human myoglobin. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L3C00301</td><td>L3C00302</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L3C00301	L3C00302
<u>Capture</u>	<u>Detection</u>						
L3C00301	L3C00302						
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.						