

SPECIFICATION SHEET

Product Name:	McAb to Myoglobin	Catalog #: L3C00302				
Description:	Monoclonal Antibody to Human Myoglobin					
Specificity:	Myoglobin					
Clone:	MY002					
Host Animal:	Mouse					
Isotype:	IgG					
Source:	Ascites					
Immunogen:	Recombinant Human Myoglobin					
Format:	Purified, Liquid					
Purification:	Protein G chromatography					
Concentration&Lot:	4.9 mg/ml	Lot: L190601				
Buffer:	10mM PBS, PH 7.4					
Preservative:	none					
Applications:	Can be used in an immunoassay for specific determination of human myoglobin. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L3C00301</td><td>L3C00302</td></tr></table>		<u>Capture</u>	<u>Detection</u>	L3C00301	L3C00302
<u>Capture</u>	<u>Detection</u>					
L3C00301	L3C00302					
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.					