

SPECIFICATION SHEET

Product Name:	McAb to CRP	Catalog #:	L3C00601
Description:	Mouse Anti Human CRP Monoclonal Antibody,		
Specificity:	Recognizes an epitope on the CRP Antigen.		
Clone:	CR001		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	Purified Human CRP Antigen		
Format:	Purified, Liquid		
Purification:	≥90% pure (SDS-PAGE).		
Concentration&Lot:	4.5mg/ml	Lot :	L180701
Buffer:	0.01 M PBS, pH 7.4 ± 0.2		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.