

## SPECIFICATION SHEET

<b>Product Name:</b>	McAb to CRP	<b>Catalog #:</b>	L3C00602
<b>Description:</b>	Mouse Anti Human CRP Monoclonal Antibody,		
<b>Specificity:</b>	Recognizes an epitope on the CRP Antigen.		
<b>Clone:</b>	CR002		
<b>Host Animal:</b>	Mouse		
<b>Isotype:</b>	IgG		
<b>Source:</b>	Ascites		
<b>Immunogen:</b>	Purified Human CRP Antigen		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	≥90% pure (SDS-PAGE).		
<b>Concentration&amp;Lot:</b>	3.7 mg/ml	<b>Lot:</b>	L190801
<b>Buffer:</b>	0.01 M PBS, pH 7.4 ± 0.2		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
<b>Storage&amp;Shelf life:</b>	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		