

SPECIFICATION SHEET

Product Name:	McAb to SAA	Catalog #:	L3C02001
Description:	Mouse Anti Serum amyloid A (SAA) Antibody,		
Specificity:	Recognizes an epitope on Serum amyloid A (SAA) Antigen.		
Clone:	SA001		
Host Animal:	Mouse		
Isotype:	IgG		
Source:	serum		
Immunogen:	Recombinant Serum amyloid A (SAA) Antigen		
Format:	Purified, Liquid		
Purification:	≥90% pure (SDS-PAGE).		
Concentration&Lot:	5 mg/ml	Lot :	L181001
Buffer:	0.01M PBS, pH 7.4 ± 0.2		
Preservative:	0.01% NaN3		
Applications:	Suitable for use in LA. Each laboratory should determine an optimum working titer for use in its particular application.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		