

## SPECIFICATION SHEET

<b>Product Name:</b>	McAb to SAA	<b>Catalog #:</b>	L3C02005
<b>Description:</b>	Mouse Anti Serum amyloid A (SAA) Antibody,		
<b>Specificity:</b>	Recognizes an epitope on Serum amyloid A (SAA) Antigen.		
<b>Clone:</b>	SA005		
<b>Host Animal:</b>	Mouse		
<b>Isotype:</b>	IgG		
<b>Source:</b>	serum		
<b>Immunogen:</b>	Recombinant Serum amyloid A (SAA) Antigen		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	≥90% pure (SDS-PAGE).		
<b>Concentration&amp;Lot:</b>	5 mg/ml	<b>Lot :</b>	L190301
<b>Buffer:</b>	0.01M PBS, pH 7.4 ± 0.2		
<b>Preservative:</b>	0.01% NaN <sub>3</sub>		
<b>Applications:</b>	Suitable for use in CLIA. Each laboratory should determine an optimum working titer for use in its particular application.		
<b>Storage&amp;Shelf life:</b>	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		