

## SPECIFICATION SHEET

**Product Name:** Recombinant FABP Protein      **Catalog #:** L4C00602

**Description:** Human Fatty Acid Binding Protein (FABP),MW:15KD

**Source:** Expressed in E. coli

**Purification:** Purified, >85% pure by SDS-PAGE

**Concentration&Lot:** 1.3mg/ml      **Lot:** L180801

**Buffer:** 20Mm Tris-HCl 140mM NaCl 2mM EDTA PH8.0

**Preservative:** None

**Applications:** Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage&Shelf life:** Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY.