

## SPECIFICATION SHEET

<b>Product Name:</b>	CMV Angiten	<b>Catalog #:</b>	L1Y00201
<b>Description:</b>	Cytomegalovirus Ag (P52+PP150) Cytomegalovirus M Antigen. (P52+PP150) WM:25KD		
<b>Source:</b>	Recombinant ,E.col		
<b>Format:</b>	Purified, Protein solution in Tris-HCl pH8.0,		
<b>Purification:</b>	Purity >85%		
<b>Concentration&amp; Lot.:</b>	2mg/ml	<b>Lot.:</b>	L180701
<b>Buffer:</b>	PBS, pH 7.4 + 0.2 , 0.5mM EDTA		
<b>Preservative:</b>	None		
<b>Applications:</b>	Immune colloidal gold immunochromatography		
<b>Storage&amp;Shelf life:</b>	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		
<b>Warnings:</b>	For research only		