

SPECIFICATION SHEET

Product Name:	SMAg-Ag	Catalog #:	L1F001
Description:	Sperm membrane antigen (SM Ag)		
Source:	Human Semen		
Format:	Ion exchange, chromatography , Purified, MW: 23-117Kd Tris buffer, pH 7.5 (5mM EDTA, 15mM mercaptoethanol, 20mM Tris).		
Concentration& Lot:	3mg/ml (Commassie Blue Staining Kit)	Lot:	L170501
Buffer:	5mM EDTA, 15mM mercaptoethanol, 20mM Tris		
Preservative:	0.2% NaN ₃		
Applications:	Suitable for use as a standard in immunoassay, Each laboratory should determine an optimum working titer for use in its particular application.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Inactivation:	Not applicable		
Warnings:	A blood sample form the tissue donor tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.