

**SPECIFICATION SHEET**

<b>Product Name:</b>	Trophoblast Antigen	<b>Catalog #:</b>	L1F004
<b>Description:</b>	Human Trophoblast Antigen (TA-Ag )		
<b>Source:</b>	Human Trophoblast		
<b>Format:</b>	Ion exchange, chromatography , Purified,Solubility is limited at neutral pH and physiological salt concentration.		
<b>Concentration&amp; Lot.:</b>	1.41mg/ml(Commassie Blue Staining Kit)	<b>Lot:</b>	L160701
<b>Buffer:</b>	20mM Tris-Hcl 150Mm NaCl ,PH7.5		
<b>Preservative:</b>	0.2% NaN3		
<b>Applications:</b>	Suitable for use as a standard in immunoassay, Each laboratory should determine an optimum working titer for use in its particular application.		
<b>Storage&amp;Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
<b>Warnings:</b>	A blood sample form the tissue donor tested negative for HIV 1, HIV 2, HCV antibodies and HBsAg and Syphilis. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.