

SPECIFICATION SHEET

Product Name:	MAb to hGH	Catalog #:	L1G00902				
Description:	Monoclonal Antibody to Human Growth Hormone (hGH)						
Clone:	GH002						
Host Animal:	Mouse						
Isotype:	IgG						
Source:	Ascites						
Immunogen:	Recombinant hGH hormone						
Format:	Purified, Liquid						
Concentration&Lot:	3.9 mg/ml (OD280nm)	Lot #:	L170201				
Buffer:	10mM PBS pH 7.4 ± 0.2						
Preservative:	None						
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00901</td><td>L1G00902</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1G00901	L1G00902
<u>Capture</u>	<u>Detection</u>						
L1G00901	L1G00902						
Storage:	Store at -20°C for at least 3 years.						

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.