

## SPECIFICATION SHEET

**Catalog #:** L2G00801                      **Lot #:** L161001

**Description:** Human Thyrogloblin (TG) Purified

**Grade:** High Purity

**Purity:** >80% by SDS-PAGE; Ion-exchange chromatography

**Source:** Thyroid tissue, surgery

**Concentration:** 6.2mg/ml

**Buffer:** 10mM PBS, pH 7.4, 0.05% NaN<sub>3</sub>.

**Biohazard:** Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.

**Storage:** Store at -20°C.