

SPECIFICATION SHEET

Product Name: IA2 **Catalog #:** L4I00501

Description: Recombinant human insulinoma associated protein 2 (IA 2)

MW: 60KD

Source: E. coli

Format: Purified, Liquid

Purification: > 85% (determined by SDS gelelectrophoresis)

Concentration& Lot.: 1.5 mg/ml **Lot.:** L180401

Buffer: 10mmPBS,1mmEDTA,pH7.4

Preservative: None

Storage&Shelf life: Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY