

SPECIFICATION SHEET

Product Name: Rabbit anti Recombinant TAU Protein **Catalog #:**L1B00603

Description: PcAb to Recombinant TAU Protein

Specificity: Reacts with Recombinant TAU Protein.

Clone: N/A

Host Animal: Rabbit

Source: Serum

Immunogen: Recombinant TAU Protein.

Format: Purified, Liquid

Purification: Protein A chromatography

Concentration&Lot: 8 mg/ml **Lot:** L190301

Buffer: 10mM PBS, pH 7.4

Preservative: none

Storage&Shelf life: Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.