

SPECIFICATION SHEET

Product Name:	Rabbit anti AD7c-NTP Protein	Catalog #: L1B00703
Description:	PcAb to Human AD7c-NTP (AD7thcloneneuronalthreadprotein) epitope peptide1 Polyclonal antibody	
Specificity:	Human AD7c-NTP epitope peptide1.	
Clone:	N/A	
Host Animal:	Rabbit	
Source:	Serum	
Immunogen:	Human AD7c-NTP epitope peptide1.	
Format:	Purified, Liquid	
Purification:	Protein A chromatography	
Concentration&Lot:	5.4 mg/ml	Lot: L190401
Buffer:	10mM PBS, pH 7.4	
Preservative:	none	
Applications:	Suitable for use in immunoassay,	
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.