

SPECIFICATION SHEET

Catalog #:	L2L009	Lot #: L190601
Description:	Procollagen Type III protein	
Source:	Human placenta	
Background:	Collagens are the major fibrillar components of most connective tissues. Procollagens are Interstitial collagens synthesised and secreted out of cells. Type III procollagen (PIIINP) is a sensitive marker of liver and lung fibrosis.	
Purification:	purified.	
Concentration:	2.36mg/ml	
Buffer:	0.5M CH ₃ COOH	
Preservative:	None	
Applications:	User optimized	
Storage:	Store at -20 °C. Avoid multiple freeze/thaw cycles.	