

SPECIFICATION SHEET

<b>Catalog #:</b>	L1V00101	<b>Lot #:</b>	L190901
<b>Description:</b>	Mouse anti-25-OH-Vitamin D3-KLH (Calcifediol-KLH) Monoclonal Antibody		
<b>Specificity:</b>	Recognizes an epitope on the 25-OH-Vitamin D3-KLH(Calcifediol-KLH)		
<b>Immunogen:</b>	25-OH-Vitamin D3-KLH (Calcifediol-KLH)		
<b>Host Animal:</b>	Mouse		
<b>Source:</b>	Serum		
<b>Format:</b>	Purified, Liquid		
<b>Concentration:</b>	1.8mg/ml		
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	ELISA		
<b>Storage:</b>	Store at -20 deg C for long term storage. Avoid repeated freeze-thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.