

SPECIFICATION SHEET

Catalog #:	LS1A01102	Lot #:	L160301				
Description:	McAb to APOE antibody Mouse anti Human APOE monoclonal Antibody,						
Specificity:	Recognizes Human APOE.						
Host Animal:	Mouse						
Immunogen:	Human APOE						
Format:	Liquid						
Buffer:	0.01M PBS, 0.15M NaCl, pH 7.4,						
Concentration:	5mg/ml						
Applications:	Suitable for use in IHC, IFA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay: <table><thead><tr><th><u>Capture</u></th><th><u>Detection</u></th></tr></thead><tbody><tr><td>LS1A01101</td><td>LS1A01102</td></tr></tbody></table>			<u>Capture</u>	<u>Detection</u>	LS1A01101	LS1A01102
<u>Capture</u>	<u>Detection</u>						
LS1A01101	LS1A01102						
Storage:	Store at -20°C. Avoid multiple freeze/thaw cycles.						

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.