

SPECIFICATION SHEET

Product Name:	Mouse anti human IgM-u chains	Catalog #:	L1100402
Description:	Mouse Antibody to human Immunoglobulin M (IgM), u chains		
Specificity:	Reacts with u chains on human IgM human immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the u region against Normal human Serum.		
Host Animal:	Mouse		
Immunogen:	Purified Human IgM(u chains)		
Format:	Ion exchange purification, Liquid		
Purification:	>90% pure (SDS-PAGE).		
Concentration:	1.3 mg/ml	Lot.:	L170201
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4		
Preservative:	None		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY