

SPECIFICATION SHEET

Product Name: Mouse IgG **Catalog #:** L11007

Description: Mouse Immunoglobulin G (IgG).

Isotype: IgG

Source: Normal mouse Serum.

Format: Liquid

Purification: >90%(SDS-page)

Concentration: 10 mg/ml **Lot:** L180901

Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4

Application: Sepecific methodologies have not been tested using this product.

Preservative: None

Storage& Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY