

SPECIFICATION SHEET

Product Name: Human IgG **Catalog #:** L1I008

Description: Human ImmunoglobulinG (IgG) Blocking Agent.

Isotype: IgG

Source: Normal Human Serum.

Immunogen: Purified Human IgG

Format: 0.01M PBS pH 7.4

Concentration: 10mg/ml **Lot:** L170801

Purification: Octoic acid purified.

Application: Sepecific methodologies have not been tested using this product.

Preservative: None

Storage& Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY