

SPECIFICATION SHEET

Product Name: Human IgM **Catalog #:** L11010

Description: Human ImmunoglobulinM (IgM)

Isotype: IgM

Source: Normal Human Serum.

Immunogen: Purified Human IgM

Concentration: 2mg/ml **Lot:** L190401

Buffer: 10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.4

Purification: Octoic acid purified.

Application: Sepecific methodologies have not been tested using this product.

Preservative: None

Storage& Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY