

SPECIFICATION SHEET

Product Name:	Goat anti Rabbit IgG (H&L)	Catalog #:	L1I01301
Description:	Goat Antibody to Rabbit Immunoglobulin G (IgG), Heavy & Light chains		
Specificity:	Reacts with heavy chains on rabbit IgG and with the light chains on all mouse immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the gamma region against Normal mouse Serum.		
Host Animal:	Goat		
Immunogen:	Purified Rabbit IgG		
Format:	Purified, Liquid		
Purification:	>90% pure (SDS-PAGE). Purified.		
Concentration:	5mg/ml	Lot.:	L181201
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2		
Preservative:	None		
Applications:	For flow cytometry and fluorescence immunohisto/cytochemistry. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY