

SPECIFICATION SHEET

Product Name:	Mouse anti human IgA	Catalog #:	L1I02101
Description:	Mouse Antibody to human Immunoglobulin A (IgA)。		
Specificity:	Reacts with heavy chains on human IgA and with the light chains on all human immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the gamma region against Normal human Serum.		
Host Animal:	Mouse		
Immunogen:	Purified human IgA		
Format:	Affinity Purified, Liquid		
Purification:	>90% pure (SDS-PAGE).		
Concentration:	1.55 mg/ml	Lot:	L161201
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2		
Preservative:	None		
Applications:	Specific methodologies have not been tested using this product.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY