

SPECIFICATION SHEET

Description: Recombinant human anti-Mullerian hormone Terminal N and C Complex

Catalog #: L2F00604

Description: Coupling Recombinant human anti-Mullerian hormone Terminal N and C Complex

Source: E. Coli **Species:** Homo sapiens

Concentration: 0.2mg/ml (Activation) **Lot #:** L190101

Purity : >90% by SDS-PAGE

Formulation : Liquid

Application: ELISA; Other applications have not been tested. The optimal dilutions should be determined by end user.

Storage buffer: 10mM CBS PH=9.0, 1mm EDTA, 5%Trehalose

Storage & Expiration: Upon receipt, aliquot and store at <-20°C
for long term. Avoid repeated freeze and thaw cycles.

Background: Müllerian inhibiting substance (MIS), also named anti-Müllerian hormone (AMH), is a tissue-specific TGF-beta superfamily growth factor. Its expression is restricted to the Sertoli cells of fetal and postnatal testis, and to the granulosa cells of postnatal ovary.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.